SOEHNLE

Body Balance Comfort F5

BEDIENUNGSANLEITUNG
OPERATING INSTRUCTIONS
MODE D' EMPLOI
ISTRUZIONI PER L'USO
GEBRUIKSAANWIJZING
INSTRUCCIONES DE MANEJO
MANUAL DE INSTRUCOES
BRUKSANVISNING
BRUGSANVISNING
KÄYTTÖOHJEET
KEZELESI UMUTATO
INSTRUKCJA OBSŁUGI
NÁVOD K POUŽITÍ
Руководство по обспуживанию
KULLAMA KILAVUZU
Οδηγίες χρήσης
NAVODILA ZA UPORABO
UPUTA ZA UPORABU

Deutsch	4
English	20
Français	36
Italiano	52
Nederlands	68
Español	84
Português	100
Dansk	116
Svenska	132
Suomi	148
Magyar	164
Polsku	180
Česky	196
по русском языку	212
Türkçe	228
Ελληνικά	244
Slovensko	260
Hrvatski	276

Deutsch

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English

Introduction	20
Important advice	20
Insert Batteries	20
Set the clock	20
Installation	21
Data input	22
Body analysis	24
Accessing your analysis values (STATE)	25
Changing personal data (SET)	26
Deleting personal data (DEL)	26
Trend curve	27
Individual recommendation	28
Factors that influence the measurement readings	31
Important notes	32
Messages	33
Technical data	33
Disposal of old electrical and electronic equipment	34
Disposal of spent batteries	34
Compliance declaration	34
Cleaning and care	34
Consumer service	34
Warranty	34

Introduction

With the Body Balance you are making a crucial contribution to your health. These electronic body analysis scales have the following functions:

- ▶ With the help of your personal data and the recommendations of health and body composition experts, Body Balance tells you your condition regarding weight, proportions of fat, water and muscle in your body and your energy consumption.
- ▶ Individual assessment of the analysis result with personal recommendations on what to do.
- ▶ Optional display of the trends of weight, proportion of fat, water and muscle in the body and energy consumption in time intervals of days, weeks or months for each standard range.
- ▶ Separate radio controlled terminal with a range of 3 metres.
- ▶ 8-person memory with automatic person recognition


Important advice

We recommend that you consult a doctor if you are overweight and are trying to reduce your weight, or if you are underweight and are trying to put on weight. Treatment and diet should only take place in consultation with your doctor.

Recommendations on exercise programmes or weight loss diets, on the basis of the values determined, should be provided by your doctor or certified health professional. Do not try to interpret the data yourself.

SOEHNLE will not be liable for damage or loss caused by the Body Balance or for third party claims.

- ▶ This product is intended for personal home use only. It is not designed for professional use in hospitals or medical facilities.
- ▶ Not suitable for people with electronic implants, such as pacemakers, etc.

Insert batteries

1. In the scales (4 x 1,5 V AA)
2. In the operator terminal (3 x 1,5 V AA)

Set the clock

3. Input hour (▼ = down, ▲ = up) and confirm (↔).
4. Input minutes (▼ = down, ▲ = up) and confirm (↔).

If you wish to change the time (Summer time /winter time) at a later point, you can switch to time mode by pressing the ▼/▲ keys at the same time.

Installation

1. Place the scales level and on a firm surface for all measurements.

Danger! Risk of slipping if surface is wet.



2. The operator terminal may optionally be held in the hand, placed on a table, remain on the scales or be hung on the wall. When fixing the magnetic wall bracket to a solid, smooth, dry and dust- and grease-free base it is sufficient for you to put the wall bracket into the oval slot on the back of the operator terminal. Then remove the tape over the adhesive patch and press the operator terminal to the fixing position horizontally at eye level.

3. With tiles or a rough surface, no guarantee can be given of the length of life of the adhesion. You should therefore screw the wall bracket in place.

To do this remove the cover from the magnetic wall bracket by inserting a fingernail or a knife point under the slot and levering off the cover.

Then drill 2 holes in the wall at a distance of 2 cm horizontally from each other. The drill size is dependent on whether you need to use wall plugs. Then screw in and press cover back in. Put the operator terminal on the magnetic bracket.

4. Option to change country-specific settings between kg/cm and lb/in or st/in. Depending on the choice of unit of weight, you can switch back to the clock mode with the  key.






The operator terminal and the scales are in radio contact.










Radio-aided computer components, HIFI components or radio-control-led toys may cause an operational error (see error messages).

Data input

To use the Soehnle Body Balance correctly, the personal data for each person need to be input into the operator terminal and you should then immediately step on the scales barefoot.



If weight is placed on the scales during input, then they immediately switch – as long as there is radio contact between the operator terminal and the scales – to the weighing and analysis mode and they cannot be used again until the scales have switched back to the clock mode (20 sec. after the weight is removed from the scales).

1. Start data input () and confirm (). If anyone's data is already stored, you must select the "NEW" display with the / keys and confirm with ().
2. Select memory (P1 ... P8) ( = down,  = up) and confirm ().
3. Input height ( = down,  = up) and confirm ().
4. Input age ( = down,  = up) and confirm ().
5. Select sex ( = male,  = female) and confirm ().


6. Input desired activity level (▼ = down, ▲ = up) and confirm (→). This level is required for determining energy consumption.

-
- 1 Scarcely any body movement (max. 2 hours/day)
-
- 2 Light, active, sitting and standing activity
-
- 3 Home and garden work, mainly standing, only occasionally seated
-
- 4 Sportspeople and people with a high level of activity
-
- 5 People involved in heavy work or high-level sport
-

7. Input the desired period for the trend display (day, week or month) (▼ = down, ▲ = up) and confirm (→). We recommend selecting a daily trend to start with because otherwise the first measurement value will not be shown in the trend display until after a week or a month.

8. Confirm "OK" (→) and, as soon as 0.0 flashes and the bar display is operating, step on the scales barefoot. If the input data are to be changed again before the first analysis, you can jump back to height with the ▼/▲ keys.









The first measurement is indispensable for storing the data for later automatic person recognition. If you don't step on the scales the message "ERROR" appears in the display and the data input must be repeated.

After you have stepped on the scales, the first body analysis will be performed as described on the following pages.

The scales have a time-out mode. It switches back to the time after about 40 seconds if the keys are not touched.

Data are retained even when batteries are changed.

Body analysis

Analysis is only possible barefoot. Personal data must have been input in advance. For meaningful results, step on the scales unclotted whenever possible and always under the same conditions.

1. Step on the scales. Your weight is then displayed and, as long as the display lasts (~5 sec.), so is a bar display. Hold still until the analysis is completed.

2. After automatic personal recognition the display once more shows the body weight, the proportions of fat (☹), water (≈), and muscles (☺) in the body and the energy consumption (**kcal**).

The normal range of the relevant person lies in the 3 middle sections of the bar display.

3. Then a personal action recommendation is displayed.

Meaning of the symbols:

- + ≈ You should drink more!
- + ☹ Muscle-building training is recommended!
- + ☺ Change of diet is recommended!
- ⚖ Your analysis figures are fine!

4. Step off the scales. The scales switch off automatically.

In rare cases, automatic assignment to a personal memory is not possible. The memories of the people who cannot be distinguished are displayed.

Since the measurement has already been taken, you can step off the scales and confirm your personal memory with the key located below it (▼ or ▲).

Then body weight, the proportions of fat, water and muscle in the body and calorie requirements are displayed.

Please note that only your body weight can be determined and no body analysis can be performed if you step on the scales in shoes or socks! The values are not entered into the trend analysis. A luxury feature of this body analysis scales product is automatic personal recognition. Recognition is achieved mainly by body weight.

In the event of changes of over ± 3 kg between individual measurements, automatic identification is no longer possible and the measurement cannot be assigned to any trend graph. This may be the case if you do not step on the scales for several days (e.g. holiday) or you weigh yourself in clothing of different weights.

Accessing your analysis values (STATE)

The latest analysis values of the people stored can be accessed on demand.


1. Start access (↵).
2. Select memoryz (P1 ... P8) with the ▼/▲ keys and confirm (↵).
3. Confirm "STATE" (↵). Then you can select the latest analysis values with the ▼/▲ keys.
4. By pressing the ↵ key again you can access the trend graphs. You can switch between the individual trend graphs with the ▼/▲ keys.

By pressing the ↵ key again you access the personal action recommendation from your last analysis.





5. By repeatedly pressing the ↵ you can choose between the different levels of the display menu.

With "EXIT" or "OFF" and confirmation with the ↵ key you can access the higher level in each case.

After 3 minutes without the keys being touched, the display switches back to the clock.

Changing personal data (SET)

If for instance:



- age is increased by one year (birthday),
- height changes (children in the growth phase)
- the activity level type is changed, then proceed as follows:

1. Start data input (←).
2. Select memory (P1 ... P8) with the ▼/▲ keys and confirm with (←).
3. "Set" with ▼/▲ keys and confirm with (←).

Then you can change the desired data as already described under "Data input".

4. After confirmation with "OK" your new data are stored. The changed data will be taken into account from the next analysis onwards.

Deleting personal data (DEL)

If data for a person – together with the trend history – are to be completely deleted, then proceed as follows:

1. Start data input (←).
2. Select memory (P1 ... P8) with the ▼/▲ keys and confirm with (←).
3. Select "DEL" with the ▼/▲ keys and confirm with (←).

A question will be asked for security purposes. Pressing the ▲ key deletes this person's entire data beyond recovery.

The memory can be reused. Quit the deletion procedure by pressing the ▼ key and the data will be retained.

Trend curve

Please weigh yourself regularly so that the trend curve will be display a precise trend.


The trend curve can display up to 59 values, subject to a corresponding number of measurements.

In case of several measurements daily, one single median value is stored. The weekly median value is determined based on 7 daily values.

The monthly median value is determined based on 30 daily values.

The value displayed after "max" and "min" is the maximum or minimum value of the curve – rounded up or down to the next whole number.


1. The calculated average range for the respective person will be displayed in the trend curve in the area between the boundary lines.
2. If a boundary line of the average value is positioned outside of the display range – in relation to the trend curve – the corresponding boundary line is displayed on the edge of the screen.
3. If only one line is displayed on the top or bottom, this means that the total curve is above or below the average range.

The display symbols of your analysis scale


* Meaning of symbols for action recommendations:

- + You should drink more fluids!
- + Muscle growth training is recommended!
- + Diet change is recommended!
- Your analysis values are ok!

Individual recommendation


The Body Balance uses your personal data and the recommendations of health experts to determine your ideal weight and body fat level.

BMI

- 1 < - 7,5
- 2 - 7,5
- 3 - 5,0
- 4 - 2,5
- 5 0,0
- 6 + 2,5
- 7 + 5,0
- 8 + 7,5
- 9 > + 7,5

Difference compared with your individual standard value

Example:

Actual body fat proportion calculated = 31.3 % of total weight.
This value is about 7 % above the calculated individual standard value.

In addition, it determines the percentage of total body water and muscle proportion. Please note that your muscles also contain body water and this is already covered by the body water analysis. Therefore the total analysis always shows a value above 100%.

Body weight

Whilst the body weight is shown in the display in kg, one of the 9 segments is shown in the bar below it.

Here the position of the active segment indicates the assessment of the body weight mentioned on the basis of the BMI figure. The BMI (Body Mass Index) is a measure for the health risk resulting from overweight or underweight.

On the basis of categorisation of the BMI, doctors distinguish between "underweight", "normal weight", "overweight" and "obese".

$$\text{BMI} = \frac{\text{body weight (kg)}}{\text{height (m)}^2}$$

Proportion of body fat

Whilst the proportion of body fat in % is shown in the display, one of the 9 segments appears in the bar below it.

Too much body fat is unhealthy and does not look nice. What is more important is that a high proportion of body fat is also usually associated with a raised proportion of fat in the blood and this increases greatly the risk of various illnesses such as for instance diabetes, heart disease, high blood pressure etc. But a very low proportion of body fat is equally unhealthy.

Apart from the fatty tissue under the skin, the body also stores important fat deposits for the protection of the internal organs and to ensure important metabolic functions.

If this essential fat reserve is compromised, this can lead to metabolic disturbance. So for instance in women with body fat figures of under 10%, menstruation may cease. In addition, the risk of osteoporosis increases. The normal proportion of body fat is dependent on age and above all on sex. In women the normal figure is about 10% higher than in men. With increasing age the body loses muscle mass and so the proportion of body fat increases.

Proportion of body water


Whilst the proportion of water in the body in % is shown in the display, one of the 9 segments appears in the bar below it.

The figure shown in this analysis reflects the "total body water" (TBW). The body of an adult is composed of about 60% water. Here there is of course a certain range, with older people showing a lower proportion of water than younger people and men showing a higher proportion than women.

The distinction between men and women is a result of the higher mass of body fat in women. Since the majority of the body's water is to be found in the fat-free mass (73% of the FFM is water) and since the water content in the fat is very low by its nature (about 10% of the fatty tissue is water), the proportion of body fat automatically decreases with an increasing proportion of water in the body.

Difference compared with your individual standard value

1	< - 12 %	very low proportion of fat
2	- 12 %	low proportion of fat
3	- 7 %	
4	- 3 %	normal proportion of fat
5	0 %	
6	+ 3 %	
7	+ 7 %	high proportion of fat
8	+ 12 %	
9	> + 12 %	very high proportion of fat


Difference compared with your individual standard value

1	< - 6 %	very low proportion of water
2	- 6 %	low proportion of water
3	- 4 %	
4	- 2 %	normal proportion of water
5	0 %	
6	+ 2 %	
7	+ 4 %	high proportion of water
8	+ 6 %	
9	> + 6 %	


Example:

Calculated actual proportion of water in the body = 51.7 % of total weight. This value is about 4 % below the calculated individual standard value.

Difference compared with your individual standard value

1	< - 6 %	very low proportion of muscle
2	- 6 %	low proportion of muscle
3	- 4 %	
4	- 2 %	normal proportion of muscle
5	0 %	
6	+ 2 %	
7	+ 4 %	high proportion of muscle
8	+ 6 %	
9	> + 6 %	

Example:

Energy consumption measured = 1,855 kcal.
This value is 5% below the calculated normal level of energy consumption.

Proportion of muscle

Whilst the proportion of muscle in % is shown in the display, one of the 9 segments appears in the bar below it.

The muscle mass is the engine of our bodies. Most of our energy conversion takes place in our muscles. The muscles take care of the maintenance of body temperature and they are the basis of our bodily performance. An exercised, muscular body is seen as aesthetically pleasing.

A check of the muscle mass is particularly important for people with weight and fat problems. Since the body needs muscles to reduce excess energy reserves or folds of fat, longterm weight reduction can only

Example:

Actual proportion of muscle calculated = 30.4 % of total body weight.
This value is about 6 % below the calculated individual standard value.

Energy consumption

While energy consumption is shown in the display in kcal, one of the 9 segments appears in the graphics panel below.

Energy consumption is calculated from the personal data – corrected by a factor for the relevant activity type - and the distribution of body mass measured.

Basic rate
+ performance rate
= energy consumption.

If you take in fewer kilocalories than your individual energy consumption, then your body will react by losing weight. Your weight will increase if you take in more calories.

Difference compared with your individual standard value

1	< - 15 %	very low energy consumption
2	- 15 %	low energy consumption
3	- 10 %	
4	- 5 %	normal energy consumption
5	0 %	
6	+ 5 %	
7	+ 10 %	high energy consumption
8	+ 15 %	
9	> + 15 %	

Factors influencing measurement readings

Analysis is based on the measurement of the body's electrical resistance. Eating and drinking habits during the course of the day and individual lifestyle affect the water balance. This is noticeable by the fluctuations in the display.

In order to ensure that the results of analysis are as accurate and consistent as possible, keep the measurement conditions constant, as only in this way will you be able to observe changes over an extended period.

Other factors can affect water balance:

- ▶ After a bath, the body fat reading may be too low and the body water reading too high.
- ▶ After a meal, readings can be higher.
- ▶ Women may experience fluctuations due to the menstrual cycle.
- ▶ Due to loss of water caused by illness or after physical activity (sport).

- ▶ After taking exercise, wait for 6 to 8 hours before carrying out the next measurement.

Varying or implausible results can occur in the case of:

- ▶ Persons with a high temperature, symptoms of oedema or osteoporosis
- ▶ Persons undergoing dialysis treatment
- ▶ Persons taking cardiovascular medicine
- ▶ Pregnant women
- ▶ Athletes who undertake more than 10 hours of intensive training per week and have a resting pulse rate of under 60/min
- ▶ Competitive athletes and body builders

Important notes

How to measure correctly?

The measurement process should be performed unclothed, if possible, and barefoot. It is recommended to perform the measurements under similar conditions and at the same time of day (e.g. in the morning, sober).

How often should measurements be performed with Soehnle body analysis scales?

Actually, that is up to you. However, we recommend always measuring at the same time under the same conditions, once a day. The daily routine, eating and drinking habits, but also your lifestyle, influence the water content of your body. When using body analysis scales, this can affect deviations in the display. The scale stores the median value.

Is it possible to recognize a change in the percentage of body mass contents on the first day?

It is possible for the body analysis scale to display deviating values. This is caused by the changing water content of your body throughout the day and the fact that your body's own values change depending on stress.

Will stored trend curves be lost if the data of one memory slot is changed?

No. Only, if the data is deleted.

Can external influences impact the measurement values?

Extreme electromagnetic influences, e.g. operating a radio transmitter in the immediate proximity of the device, can influence the measurement results.


After removal of the influence, the product will be fully functional again.

However, it might be necessary to switch the device off and back on.

How are trend curves created?

Daily trend curve (day) means:

The daily median values of the last 59 days will be displayed. The last value will always be added to the right side.

Weekly curves (week) means:

A weekly median value will be determined based on 7 daily median values. Therefore, the first weekly median value can be displayed only after one week*.

The last 59 weekly median values will be displayed.

Monthly curve (monthly) means:

A monthly median value will be determined based on 30 daily median values. Therefore, the first monthly median value can be displayed only after one month*. Up to 59 monthly median values can be displayed.

* We therefore recommend to select a daily trend initially.

Messages

1. Batteries in the operator unit (LOW BAT-T) or the base (LOW BAT-B) low.

2. Overload: From 150 kg.

3. Implausible analysis result – check individual programming.

4. + 5.

Bad foot-contact:
clean scale or feet or check for sufficient skin moisture.

6. If skin is too dry – moisten feet or carry out measurement after taking a bath or shower. Interference with radio contact between scales and operator terminal.

The range of the radio signal and correct operation are dependent on several factors:


Battery condition, location and position, screening objects, other devices on the same frequency band.

Under the environment with electrostatic discharge, the 'Body Balance Comfort F5' could not maintain communication link and require the user to take off the battery and insert it again.

Under the environment with radio interference the 'Body Balance Comfort F5' could not maintain communication link and resume to normal when the interference stops.

Technical data

- ▶ Weighing capacity x graduation = max. 150kg x 100g
- ▶ Display resolution: body fat: 0.1%
- ▶ Display resolution: body water: 0.1%
- ▶ Display resolution: muscle mass: 0.1%
- ▶ Display resolution: energy consumption: 1 kcal
- ▶ height: 100–250 cm
- ▶ age: 10–99 years
- ▶ 8-person individual memory locations
- ▶ Size of LCD screen
- ▶ measuring circuit: 0.035 mA RMS, 50 kHz
- ▶ Radio frequency: 433 MHz, range: 3m
- ▶ Batteries: 4 x 1,5V AA for the scales
3 x 1,5 V AA for the operator terminal
Batteries supplied with unit.

Disposal of old electrical and electronic equipment

This symbol on the product or on its packaging indicates that this product shall not be treated as household waste.

Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. For more detailed information about recycling contact your local city office, your household waste disposal service or the shop where you purchased the product.

Disposal of spent batteries

Batteries must not be disposed of as normal household waste. Note that you are legally obliged to dispose of used batteries in a correct way. You can return spent batteries either to public collection points in your town or to any outlet selling batteries of the same kind.

Pb = Battery contains lead

Cd = Battery contains cadmium

Hg = Battery contains mercury

Compliance declaration

This device is screened in accordance with the applicable EC Directive 1999/5/EC.

Note: the display value may be affected by extreme electromagnetic influences, e.g. when a radio is operated in the immediate vicinity of the device.

The product can be used for its intended purpose again when the interference disappears (a reset may be required).

Cleaning and care

Only with a slightly damp cloth.

Do not use any solvents or abrasive substances and do not immerse the scales in water.

Consumer service

If you have questions or comments, the following points of contact are available:

Ph: (08 00) 5 34 34 34

Warranty

SOEHNLE guarantees that all defects due to materials or manufacturing faults will be remedied by replacement or repair, free of charge, for a period of 3 years from the date of purchase. Please keep your purchase receipt and the guarantee card in a safe place.

If you have any complaints, please return the scale to your dealer with the guarantee card and receipt.

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LEIFHEIT AG