UV Bracelet

Art.no 31-4819

Please read the entire instruction manual before using the product and save it for future reference. We reserve the right for any errors in text or images and any necessary changes made to technical data. If you have any questions regarding technical problems please contact our customer services.

Safety

- This product contains no sunblock and offers no protection from sunlight.
- The use of this product for limited sunbathing in combination with sun cream and clothing can reduce the risk of skin damage, but not prevent it entirely.
- Unprotected skin which has been repeatedly exposed to sunlight can age prematurely. There is still a risk of contracting skin cancer and other ailments even if you have not been sunburned.
- Use just one bracelet per person per day and discard the bracelet once it has been used.
- The bracelet only works with SPF15 or higher.

Use

1. Open the packet, take out a bracelet and then close the packet immediately.
2. Fasten the bracelet around your wrist with the printed side facing outwards.
3. Apply sunscreen to all exposed body parts and over the bracelet.
4. When the bracelet is exposed to sunlight, it will turn purple to indicate that it has been activated.
5. When the bracelet turns pink it means that it is time to apply more sunscreen to the bracelet and your skin.
6. When the bracelet turns yellow you should move into the shade or put on some clothes.

(For color reference check the packaging.)

When a partly used bracelet is exposed to sunlight again, its indication level will resume from the time it was removed from the sunlight.

Disposal

Follow local ordinances when disposing of this product. If you are unsure how to proceed, contact your local authority.

Tilda's AB
Mejselgatan 12
S- 235 32 Vellinge

Hereby declares that this product conforms to *Direktiv 93/42/EEC* valid until 2012-12-31

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