Massage Unit

ENGLISH

Art.no: 18-4235 Model: JL-3003-UK
34-3012 JL-3003

Please read the entire instruction manual before using and save it for future reference. We apologise for any text or photo errors and any changes of technical data. If you have any questions concerning technical problems please contact our Customer Service Department.

- **Deep penetrating massage. Provides effective relief of tension and muscle pain.**
- **Variable speed control and ergonomic handle.**
- **2 accessories for various types of massage included.**

Safety

Caution: Do not open the casing. Certain components inside the casing carry dangerous currents. Contact with these can give electric shock.

1. This product must not be used by anyone (including children) suffering from physical or mental impairment. This product should not be used by anyone who does not have the proper experience or knowledge of the product unless they have been instructed in its use by someone who takes responsibility for their safety.
2. Never let children play with the product.
3. Always pull the plug after use. Never leave the product unattended when plugged into the mains.
4. Never use the massager in the bathroom or in damp environments.
5. Do not place the massager where it could fall into water or any other type of fluid.
6. An electric appliance which is connected to a wall socket is still live even when not in use.
7. Never touch the massager with wet hands when it is connected to the mains.
8. If the massager is connected to the mains and falls into water: **PULL OUT THE PLUG** before touching the massager.
9. Pull out the plug before cleaning.
10. Never use the massager if the plug, mains lead or any other part is damaged.
11. Only use the supplied charger. The product should not be taken apart or modified.
12. For indoor use only. Do not use the product in environments where there are flammable gases or vapours present such as propellants, oxygen cylinders, etc.
13. Make sure that no foreign objects come into contact with any openings on the product as this can lead to electric shock.
14. The power lead must be changed if damaged. This is to prevent the risk of electric shock or fire and should only be carried out by the manufacturer, qualified service facility or qualified technician.
15. The product should only be repaired by qualified service personnel.

WARNING: Do not use the product in the proximity of the bath, shower, washbasin or other containers that contain liquid.

Operation

Caution:

Seek a doctor's advice before using the device if you are pregnant.

Do not massage more than 5 minutes on the same area or more than 30 minutes at a time.

1. Make sure the power switch is in the **OFF** position. Plug in the massager.
2. Set the switch to the **ON** position. Adjust the speed with the dial from mild (**MIN**) to intensive (**MAX**).
3. When finished, turn the power switch to the **OFF** position and unplug the massager.

Tip: Drink a few glasses of water before and after massaging since this assists in extracting waste from the body.

Care

Clean the massager every now and then with a lightly moistened cloth. Use a mild detergent and avoid solvents and strong, abrasive detergents.

Recycling

Follow local ordinances when disposing of this product. If you are unsure of how to dispose of this product, please contact your municipality.

Specifications

Model: JL-3003, JL-3003-UK
Rated Voltage: 220/240 V AC, 50/60 Hz
Power: 35 W
CE approved.

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